

Synonym


CUB domain-containing protein 1,CDCP1,CD318,Membrane glycoprotein gp140,Subtractive immunization M plus HEp3-associated 135 kDa protein,SIMA135,Transmembrane and associated with src kinases,TRASK,UNQ2486,PRO5773

Source

Biotinylated Human CDCP1 Protein, His,Avitag(CD1-H82E4) is expressed from human 293 cells (HEK293). It contains AA Phe 30 - Thr 667 (Accession # Q9H5V8-1).

Predicted N-terminus: Phe 30

Molecular Characterization



This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (AvitagTM)

The protein has a calculated MW of 75.7 kDa. The protein migrates as 90-110 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using AvitagTM technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20 $^{\circ}$ C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Biotinylated Human CDCP1 Protein, His, Avitag on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

SEC-MALS



The purity of Biotinylated Human CDCP1 Protein, His, Avitag (Cat. No. CD1-H82E4) is more than 85% and the molecular weight of this protein is around 90-110 kDa verified by SEC-MALS.

<u>Report</u>

Biotinylated Human CDCP1 Protein, His,Avitag™ (MALS verified)

Catalog # CD1-H82E4



Background

CUB domain-containing protein 1, also know as CDCP1. CDCP1 is a cell surface glycoprotein that has been recognized both as a tumour marker and as a potential target to disrupt progression of cancer. CDCP1 has been shown to modulate cellsubstratum adhesion and motility in colon cancer cell lines.

Clinical and Translational Updates

Please contact us via <u>TechSupport@acrobiosystems.com</u> if you have any question on this product.