



Source

Biotinylated Mouse H-2Kd&B2M&InsB (LYLVCGERL) Complex Protein(H2B-M82E3) is expressed from human 293 cells (HEK293). It contains AA Gly 22 - Thr 305 (H-2Kd) & Ile 21 - Met 119 (B2M) & LYLVCGERL peptide (Accession # [P01902](#) (H-2Kd) & [P01887](#) (B2M) & LYLVCGERL). Predicted N-terminus: Gly 22 & Ile 21

Molecular Characterization

Biotinylated Mouse H-2Kd&B2M&InsB (LYLVCGERL) Complex Protein is produced by co-expression of H-2Kd and B2M loaded with InsB peptide.

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

The protein has a calculated MW of 40.0 kDa and 15.6 kDa. The protein migrates as 50-55 kDa and 17 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Biotinylated Mouse H-2Kd&B2M&InsB (LYLVCGERL) Complex Protein on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS



The purity of Biotinylated Mouse H-2Kd&B2M&InsB (LYLVCGERL) Complex Protein (Cat. No. H2B-M82E3) is more than 90% and the molecular weight of this protein is around 55-75 kDa verified by SEC-MALS. [Report](#)

Discounts, Gifts,
and more!





Clinical and Translational Updates

**Discounts, Gifts,
and more!**

