

Synonym

HLA-A*0301 & B2M & KRASG12C (VVVGACGVGK)

Source

Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein(HLC-H52H9) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Thr 305 (HLA-A*03:01) & Ile 21 - Met 119 (B2M) & VVVGACGVGK peptide (Accession # [P04439](#) (HLA-A*03:01) & [P61769](#) (B2M) & VVVGACGVGK).

Predicted N-terminus: Gly 25 & Ile 21

Molecular Characterization

Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein is assembled by biotinylated monomer (HLC-H82E3) and streptavidin.

Biotinylated Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Complex Protein is produced by co-expression of HLA and B2M loaded with KRASG12C peptide. Biotinylated Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Complex Protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

The protein has a calculated MW of 35.9 kDa and 11.7 kDa. The protein migrates as 40-43 kDa, 13 kDa and 10 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE


Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA

Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein ELISA
0.1 µg of Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein per well




Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein ELISA
0.1 µg of Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein per well



Immobilized Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein (Cat. No. HLC-H52H9) at 1 µg/mL (100 µL/well) can bind Anti-HLA class I Antibody, Human IgG1 (W6/32) with a linear range of 0.1-2 ng/mL (QC tested).

Immobilized Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein (Cat. No. HLC-H52H9) at 1 µg/mL (100 µL/well) can bind Anti-B2M Antibody, Human IgG1 with a linear range of 0.1-2 ng/mL (Routinely tested).

Bioactivity-SPR



Human HLA-A*03:01&B2M&KRASG12C (VVVGACGVGK) Tetramer Protein (Cat. No. HLC-H52H9) captured on CM5 Chip via Anti-B2M antibody can bind Anti-HLA class I Antibody, Human IgG1 (W6/32) with an affinity constant of 2.31 pM as determined in a SPR assay (Biacore 8K) (Routinely tested).

Background

The Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) oncogene plays a critical role in the initiation and maintenance of pancreatic tumors and its signaling network represents a major target for therapeutic intervention. The Human HLA-A*0301 KRASG12C (VVVGACGVGK) complex protein is a complex of HLA-A*0301 of the MHC Class I, B2M, and VVVGACGVGK peptide of the KRASG12C.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.