

PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein

Catalog # HLV-HP2H7



BIOSYSTEMS
Acro

Synonym

HLA-A*1101 & B2M & KRASG12V (VVVGAVGVGK)

Source

PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein(HLV-HP2H7) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Thr 305 (HLA-A*11:01) & Ile 21 - Met 119 (B2M) & VVVGAVGVGK peptide (Accession # [Q5S3G3-1](#) (HLA-A*11:01) & [P61769](#) (B2M) & VVVGAVGVGK).

Predicted N-terminus: Gly 25 & Ile 21 & Val

Molecular Characterization

PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein is assembled by biotinylated monomer (HLV-H82E7) and PE-labeled streptavidin.

Biotinylated Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Complex Protein is produced by co-expression of HLA and B2M loaded with KRASG12V peptide. Biotinylated Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Complex Protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

Conjugate

PE

Excitation Wavelength: 488 nm / 561 nm

Emission Wavelength: 575 nm

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, 1% BSA, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Please protect from light and avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


Bioactivity-ELISA

PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein ELISA
0.1 µg of PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein per well



Immobilized PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein (Cat. No. HLV-HP2H7) at 1 µg/mL (100

PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein ELISA
0.1 µg of PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein per well



Immobilized PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein (Cat. No. HLV-HP2H7) at 1 µg/mL (100

Discounts, Gifts,
and more!



PE-Labeled Human HLA-A*11:01&B2M&KRASG12V (VVVGAVGVGK) Tetramer Protein

Catalog # HLV-HP2H7



BIOSYSTEMS
Acro

μL/well) can bind Anti-HLA class I Antibody, Human IgG1 (W6/32) with a linear range of 0.12-1 ng/mL (QC tested).

μL/well) can bind Anti-B2M Antibody, Human IgG1 with a linear range of 0.12-1 ng/mL (Routinely tested).

Background

The Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) oncogene plays a critical role in the initiation and maintenance of pancreatic tumors and its signaling network represents a major target for therapeutic intervention. The Human HLA-A*1101 KRASG12V (VVVGAVGVGK) complex protein is a complex of HLA-A*1101 of the MHC Class I, B2M, and VVVGAVGVGK peptide of the KRASG12V.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

2/4/2024