

Product Details

Laminin 411 is a key component of the basement membrane and is involved in the regulation of cell differentiation. Laminin 411 fragment consists of the Cterminal end of the alpha 4, beta 1, and gamma 1 chains.

Key parameter

Purity (SDS PAGE) > 90%

Mycoplasma Test Negative

Sterility Test Negative

Endotoxin Test < 0.1 EU per μg

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 24 months in lyophilized state;
- -70°C for 12 months under sterile conditions after reconstitution;
- 2-8°C for 3 months under sterile conditions after reconstitution.


Bioactivity-ELISA

Human Laminin 411 Protein, premium grade ELISA 0.5 μg of Human ITGA6&ITGB1 Heterodimer Protein, His Tag&Tag Free per well



Immobilized Human ITGA6&ITGB1 Heterodimer Protein, His Tag&Tag Free (Cat. No. IT1-H52W7) at 5 μ g/mL (100 μ L/well) can bind Human Laminin

Human Laminin 411 Protein, premium grade ELISA 0.5 μg of Human ITGA7X2B&ITGB1 Heterodimer Protein, His Tag&Tag Free per well



Immobilized Human ITGA7X2B&ITGB1 Heterodimer Protein, His Tag&Tag Free (Cat. No. IT1-H52W8) at 5 μg/mL (100 μL/well) can bind Human







411 Protein, premium grade (Cat. No. LAT-H5263) with a linear range of 0.1-2 ng/mL (QC tested).

Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) with a linear range of 0.01-0.625 μg/mL (Routinely tested).


Bioactivity-SPR




Human Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) immobilized on CM5 Chip can bind Human ITGA6&ITGB1 Heterodimer Protein, His Tag&Tag Free (Cat. No. IT1-H52W7) with an affinity constant of 29.3 nM as determined in a SPR assay (Biacore 8K) (QC tested).




Human MCAM, His Tag (Cat. No. CD6-H5229) immobilized on CM5 Chip can bind Human Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) with an affinity constant of 0.115 μ M as determined in a SPR assay (Biacore 8K) (Routinely tested).



Human Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) captured on Protein A Chip can bind Human MCAM, His Tag (Cat. No. CD6-H5229) with an affinity constant of 107 nM as determined in a SPR assay (Biacore 8K) (Routinely tested).



Human Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) immobilized on CM5 Chip can bind Human ITGA7X2B&ITGB1 Heterodimer Protein, His Tag&Tag Free (Cat. No. IT1-H52W8) with an affinity constant of 17.4 nM as determined in a SPR assay (Biacore 8K) (Routinely tested).




Human MCAM, His Tag (Cat. No. CD6-H5229) captured on CM5 chip via anti-His antibody can bind Human Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) with an affinity constant of 27.9 nM as determined in a SPR assay (Biacore 8K) (Routinely tested).

Human Laminin 411 Protein, premium grade

Catalog # LAT-H5263



Bioactivity-BLI



Loaded Human Laminin 411 Protein, premium grade (Cat. No. LAT-H5263) on Protein A Biosensor, can bind Human ITGA6&ITGB1 Heterodimer Protein, His Tag&Tag Free (Cat. No. IT1-H52W7) with an affinity constant of 26.9 nM as determined in BLI assay (ForteBio Octet Red96e) (QC tested).

Clinical and Translational Updates

